

LIVINGSTONE INSULIN SYRINGES

For accurate and comfortable administration of prescribed insulin through subcutaneous injection for medication process

- Sterile
- White plunger
- With needle
- Used with U-100 insulin only
- Individually packed

Unit: 100 Pieces/Box

Item No.	Capacity	Needle Size
DSL0003ML29G	0.3mL	29G x 0.5"/12.7mm
DSL0003ML31G	0.3mL	31G x 0.32"/8mm
DSL0005ML27G	0.5mL	27G x 0.5"/13mm
DSL0005ML29G	0.5mL	29G x 0.5"/12.7mm
DSL0005ML30G	0.5mL	30G x 0.32"/8mm
DSL0005ML31G	0.5mL	31G x 0.32"/8mm
DSL001ML27G	1mL	27G x 0.5"/12.7mm
DSL001ML29G	1mL	29G x 0.5"/12.7mm
DSL001ML30G	1mL	30G x 0.32"/8mm
DSL001ML31G	1mL	31G x 0.32"/8mm

Unit: Each

Item No.	Capacity	Needle Size
DSL0003ML31GL	0.3mL	31G x 0.32"/8mm
DSL0005ML30GL	0.5mL	30G x 0.32"/8mm
DSL001ML31GL	1mL	31G x 0.32"/8mm



Item No.: **DSL0003ML31G**



Item No.: **DSL0005ML31G**



Item No.: **DSL001ML29G**

BUY FROM OUR WEBSTORE WWW.LIVINGSTONE.COM.AU
AND GET THE LATEST RANGE AND BEST PRICES



Follow us on **INSTAGRAM**
[@livingstoneinternational](https://www.instagram.com/livingstoneinternational)

