



INTRAVENOUS (IV) CATHETERS

Sample blood and administer fluids with ease

- Use of hydrophobic membrane filter eliminates blood leakage
- Can be connected to syringe by removing filter cap to expose Luer taper end
- Pyrogen free
- Sterile - Sterilised by ETO (Ethylene Oxide Processing)
- Transparent catheter hub provides confirmation of successful catheter placement when withdrawing the needle
- Colour coded for easier identification of catheter size

Unit: **50 Pieces/Box**



14G x 2 Inches
Item No.: **CIV1451N**



16G x 2 Inches
Item No.: **CIV1651N**



18G x 2 Inches
Item No.: **CIV1851N**



20G x 1 Inch
Item No.: **CIV2025N**



22G x 1 Inch
Item No.: **CIV2225N**

BUY FROM OUR WEBSTORE WWW.LIVINGSTONE.COM.AU
AND GET THE LATEST RANGE AND BEST PRICES



Follow us on **INSTAGRAM**
[@livingstoneinternational](https://www.instagram.com/livingstoneinternational)

