

## Reliable Control for Safer Blood Collection

Designed for dependable performance in clinical settings, the **Livingstone Tourniquet** assists with controlled venous restriction during blood collection, infusion and transfusion procedures. Its elastic fabric construction provides a secure yet comfortable fit, while the buckle release design allows quick application and removal. Suitable for reusable single patient use.



### Key Features:

- **Elastic Fabric Construction:** Provides flexibility and comfortable application around the limb.
- **Buckle Release Design:** Allows quick fastening and easy removal during procedures.
- **Hook Fastener Patches:** Helps maintain a secure fit during clinical use.
- **Latex Free Material:** Reduces the risk of latex-related sensitivities.
- **Reusable Format:** Supports cost-effective use for individual patients.
- **Soft Plush Surface:** Enhances patient comfort during application.

## Reliable Control for Safer Blood Collection

### Intended Use:

- Assists in slowing venous blood flow during blood collection procedures
- Suitable for infusion and transfusion preparation
- Designed for use in hospitals, clinics, pathology centres and healthcare facilities
- Intended for reusable single patient use

Item No.	Description	Unit
1 <b>TOURBUCKR</b>	Livingstone Tourniquet, 43 x 2.5cm, Elastic Fabric, Buckle Release, Latex Free, Reusable, Blue	25 Pieces/Box
2 <b>TOURBUCKC</b>	Livingstone Tourniquet, 43 x 2.5cm, Elastic Fabric, Buckle Release, Latex Free, Reusable, Blue	1 Piece/Blister Card
<b>TOURBUCK</b>	Livingstone Tourniquet, 43 x 2.5cm, Elastic Fabric, Buckle Release, Printed Single Patient Use, Latex Free, Blue	25 Pieces/Box
<b>TOURBUCKL</b>	Livingstone Tourniquet, 43 x 2.5cm, Elastic Fabric, Buckle Release, Printed Single Patient Use, Latex Free, Blue	Each

Built for busy healthcare environments, the **Livingstone Tourniquet** supports efficient patient preparation while helping maintain clinician workflow. The practical fastening system and durable elastic weave make it a dependable addition to everyday procedural settings where consistency and ease of handling matter.

### Product Specification:

- **Material:** Elastic Woven Fabric
- **Design:** Buckle Release with Adjustable Fastening
- **Dimensions:** 43 x 2.5cm
- **Reusability:** Reusable Single Patient Use
- **Pack Type:** Available in Individual and Box
- **Colour:** Blue
- **Latex Free:** Yes

### Care:

- For single patient reuse only
- Inspect before each use for signs of wear or damage
- Do not apply excessively tight pressure
- Clean according to facility infection control protocols
- Store in a cool, dry place
- Discontinue use if material deterioration occurs

Buy from our webstore and get the latest range and best prices  
[www.livingstone.com.au](http://www.livingstone.com.au)