
BLOOD COLLECTION TUBES

LIVINGSTONE


Ensuring Secure Collection and Dependable Sample Integrity

Livingstone Blood Collection Tubes are high-quality laboratory tubes designed for the efficient collection, preservation, and transport of blood samples, supporting accurate diagnostic testing across a wide range of clinical and laboratory applications.


1



2



3



4





Key Features:

- **Comprehensive Tube Range:** Includes multiple tube types with specific additives such as EDTA, Sodium Citrate, Fluoride/Oxalate, and Serum Separators.
- **Colour-Coded Cap System:** Standardised cap colours allow quick identification of tube type, helping streamline workflow and reduce handling errors.
- **Secure Closure Mechanisms:** Features non-spill caps, Hemogard closures, and pull caps to enhance safety and minimise leakage risks.
- **Accurate Additive Ratios:** Manufactured to ensure correct blood-to-additive proportions for reliable and consistent test outcomes.
- **Material Versatility:** Available in both plastic and glass tubes to accommodate different laboratory preferences and procedures.
- **Laboratory Compatibility:** Designed for seamless use with standard blood collection systems and diagnostic analysers.

livingstone.com.au sales@livingstone.com.au 1300 548 289 *(LIV BUY)*

BLOOD COLLECTION TUBES

LIVINGSTONE

Ensuring Secure Collection and Dependable Sample Integrity

Intended Use:

- For collection and containment of venous blood specimens
- Suitable for hematology, coagulation, chemistry, and serology testing
- Designed to preserve sample integrity during handling and transport
- Intended for use by trained healthcare professionals in clinical and laboratory environments

Item No.	Description	Unit
BD36903	Livingstone Blood Collection Glass Tube, 4mL, 13 x 75mm, Red Hemogard Closure, Serum Tubes, Plain	100 Pieces/Pack
BD361017	Livingstone Blood Collection Glass Tube, 5mL, 13 x 75mm, Pink Hemogard Closure, Aprotinine K3 EDTA Tube	100 Pieces/Pack
BD365975	Livingstone Blood Collection Plastic Tube, 0.5mL, 10 x 45mm, Pink Closure, K2 EDTA Additive	50 Pieces/Rack
GRN454034	Livingstone Blood Collection Plastic Tube, 1mL, 13 x 75mm, Lavender Cap, K3 EDTA Tube	50 Pieces/Rack
BD366664	Livingstone Blood Collection Plastic Tube, 2mL, 13 x 75mm, Green Hemogard Closure, Lithium Heparin Tube	100 Pieces/Rack
1 BCTCIT	Livingstone Blood Collection Plastic Tube, 2.7mL, 13 x 75mm, Light Blue Non-Spill Cap, 3.2% Sodium Citrate	100 Pieces/Pack
EDTA454020	Livingstone Blood Collection Plastic Tube, 3mL, 13 x 75mm, Lavender Cap, K2 EDTA Tube	50 Pieces/Rack
2 BCTEDTA	Livingstone Blood Collection Plastic Tube, 4mL, 13 x 75mm, Purple Lavender Non-Spill Cap, K3 EDTA	100 Pieces/Pack
BD367839	Livingstone Blood Collection Plastic Tube, 4mL, 13 x 75mm, Lavender Hemogard Closure, K2 EDTA	100 Pieces/Pack
BD367883	Livingstone Blood Collection Plastic Tube, 4mL, 13 x 75mm, Green Hemogard Closure, Lithium Heparin Tube	100 Pieces/Rack
BCTFLOX	Livingstone Blood Collection Plastic Tube, 5mL, 13 x 100mm, Grey Non-Spill Cap, Fluoride/Oxalate	100 Pieces/Pack
BCTSST5	Livingstone Blood Collection Plastic Tube, 5mL, 13 x 100mm, Gold Non-Spill Cap, SST, Clot Activator, Gel Separator & Coagulant	100 Pieces/Pack
BCTEDTA6	Livingstone Blood Collection Plastic Tube, 6mL, 13 x 100mm, Royal Blue Non-Spill Cap, K2 EDTA	100 Pieces/Pack
3 BD367885	Livingstone Blood Collection Plastic Tube, 6mL, 13 x 100mm, Green Hemogard Closure, Lithium Heparin Tube	100 Pieces/Rack
BD367876	Livingstone Blood Collection Plastic Tube, 6mL, 13 x 100mm, Green Hemogard Closure, Sodium Heparin Tube	100 Pieces/Rack
GBO456001	Livingstone Blood Collection Plastic Tube, 6mL, 13 x 100mm, White Pull Cap, No Additive, Non-Ridged	50 Pieces/Rack
BD367377	Livingstone Blood Collection Plastic Tube, 8mL, Light Green Hemogard, Plasma Separator Tube, Lithium Heparin and Gel, Sterile, PET	100 Pcs/Rack

BLOOD COLLECTION TUBES

LIVINGSTONE

Ensuring Secure Collection and Dependable Sample Integrity

Item No.	Description	Unit
4 **BCTSST**	Livingstone Blood Collection Plastic Tube, 8.5mL, 16 x 100mm, Gold Non-Spill Cap, SST, Clot Activator, Gel Separator & Coagulant	100 Pieces/Pack
BD366645	Livingstone Blood Collection Plastic Tube, 8.5mL, 16 x 100mm, Yellow Conventional Closure, ACD Solution A Tube	100 Pieces/Rack
BCTEDTA9	Livingstone Blood Collection Plastic Tube, 9mL, 16 x 100mm, Purple Lavender Non-Spill Cap, K3 EDTA	50 Pieces/Pack
GBO455001	Livingstone Blood Collection Plastic Tube, 9mL, 16 x 100mm, White Pull Cap, No Additive	50 Pieces/Pack
BCTSER	Livingstone Blood Collection Plastic Tube, 10mL, 16 x 100mm, Red Non-Spill Cap, Serum Tube	100 Pieces/Pack
BCTEDTAK2	Livingstone Blood Collection Plastic Tube, 10mL, 16 x 100mm, Purple Lavender Non-Spill Cap, K2 EDTA	100 Pieces/Pack

Livingstone Blood Collection Tubes are designed to support smooth and efficient laboratory workflows from sample collection through to analysis. Their standardised format helps reduce handling variability, allowing healthcare professionals to maintain consistency across high-volume testing environments. By supporting reliable specimen quality and streamlined processing, they contribute to improved operational efficiency in both clinical and laboratory settings.

Product Specifications:

- **Tube Material:** Medical-Grade Plastic or Glass
- **Tube Sizes:** Available in 1mL, 2mL, 2.7mL, 3mL, 4mL, 5mL, 6mL, 8mL, 8.5mL, 9mL, 10mL
- **Additives:** Anticoagulant Tubes, Aprotinine, Clot Activator, Fluoride/Oxalate, Gel Separator, K2 EDTA, K3 EDTA, Lithium Heparin, No Additive, Plain Tubes, Plasma Separator, Serum Tubes, Sodium Citrate (3.2%), Sodium Heparin, SST (Serum Separator Tubes)
- **Closure Types:** Non-Spill Caps, Hemogard Closures, Pull Caps
- **Cap Colors:** Light Blue, Purple/Lavender, Red, Pink, Green, Light Green, Grey, Gold, Royal Blue, White

Care:

- Store in a clean, dry place
- Handle carefully to avoid breakage, especially glass tubes
- Do not use if packaging is damaged or compromised
- Keep sealed until ready for use
- Follow standard laboratory biosafety procedures
- Dispose of used tubes in biohazard waste containers
- Avoid exposure to extreme temperatures
- Check expiry date before use

Buy from our webstore and get the latest range and best prices
www.livingstone.com.au